
INSTRUCTIONS FOR USE
NAUDOJIMO INSTRUKCIJOS
LIETOŠANAS INSTRUKCIJA
KASUTUSJUHEND

DOLORCLAST[®]

FOCUSED SHOCK WAVES

CE₀₁₂₄

EMS⁺

ENGLISH.....	3
LIETUVIŲ	31
LATVIEŠU	59
EESTI.....	87

Images are non contractual and indicative. These may evolve without notice.
Paveikslėliai nėra sutarties objektas, o tik orientaciniai. Jie gali būti keičiami be išankstinio įspėjimo.
Attēli ir ārpustīgumiski. Tie var tikt attīstīti bez brīdinājuma.
Pildid on mittelepingulised ja näitlikud. Neid võidakse ilma ette teatamata muuta.

CONTENTS

ENGLISH

LIETUVIŲ

LATVIEŠU

EESTI

1. Please read this first	4	5. Getting started	16
1.1. About these Instructions for Use	4	5.1. Starting up the console	16
1.2. Description of the DolorClast® Focused Shock Waves	4	5.2. Interface	16
1.3. Operating principle of the handpiece	4	5.2.1. Configure your background	16
1.4. Intended use	4	5.2.2. General modes	16
1.5. Indications for use	5	5.2.3. Special modes: Analgesic and Ramp up	16
1.6. Intended clinical benefits	5	5.2.4. VAS assessment and Tutorials	17
1.7. Intended user and environment	5	6. Treatment	18
1.8. Intended patient population	5	6.1. Prerequisites	18
1.9. Contraindications	5	6.2. Preparing the patient and the handpiece	18
1.10.Side-effects	5	6.3. Preset protocols and treatment settings	19
1.11.Potential complications	6	6.3.1. Manual settings: Intensity / Energy, Frequency, Impulses	19
1.12.Warnings and precautions	6	6.3.2. Indication screen	19
1.13.Indications	7	6.3.3. Factory settings: Indications and preset protocols	20
2. Warning	8	6.3.4. Creating or Saving your custom protocol	21
3. Components	8	6.4. Performing the treatment	21
3.1. Figure 1: Unpacking	8	6.4.1. Handling the handpiece	21
3.2. Figure 2: Handpiece and gel pad kit	9	6.4.2. Starting the treatment	22
3.3. Figure 3: Console	9	6.4.3. Assessing patient pain	22
3.4. Touchscreen interface introduction – Main HMI screens	10	6.4.4. Pausing / Ending the treatment	22
3.4.1. The Ready-to-Work screens: Treatment, Indication, Custom	10	7. Reprocessing and Maintenance ..	23
3.4.2. SETTINGS Screens: Console - Treatment ..	12	7.1. Cleaning and disinfecting	23
3.4.3. Factory default settings	13	7.2. Maintenance	23
4. Installation	14	8. Product disposal	24
4.1. Installing the console	14	9. Troubleshooting - Warranty - Support	24
4.1.1. Installing the console on the cart	14	9.1. Troubleshooting table	24
4.1.2. Installing the handpiece holder onto the device	14	9.2. Error codes	25
4.2. Connecting the console to the equipotential conductor	15	9.3. Warranty	25
4.3. Connecting the handpiece to the console	15	9.4. EMS Technical Support	25
4.4. Connecting the power cable	15	10.Symbols / Glossary	26
		10.1.Symbols	26
		10.2.Glossary	26
		11.Technical description	27
		12.Electromagnetic compatibility ...	28

To make these Instructions for Use simpler to read, the brand names are written as follows:

DolorClast® Focused Shock Waves is your product

DOLORCLAST® FOCUS handpiece, as the emission source, is the handpiece

Focused Shock Waves Console is the console

Compact Cart DolorClast® is the cart (optional)

1. PLEASE READ THIS FIRST

1.1. About these Instructions for Use

Thank you for purchasing a DolorClast® Focused Shock Waves.

We are more than happy to answer questions you may have, and we welcome your feedback and suggestions. We are also on-hand to provide support in the event of any technical issues.

Please contact your EMS authorized service center or your distributor directly.

We wish you lots of success!

Please note that the English version of this manual is the source from which all translations are derived. In case of any discrepancy, the binding version is the English text.

Before using the DolorClast® Focused Shock Waves, please read carefully and follow the recommendations set out in the Instructions for Use originally delivered with your device and system.

Please pay close attention to the safety precautions.

Any serious incident that occurs in relation to the product should be reported to the manufacturer and the competent authority.

For technical issues or complaints related to the use of this product, please contact your local customer service representative.

To prevent injury to people and damage to property, please follow the corresponding instructions, marked as indicated below.

WARNING

Risk of severe patient or user injuries

CAUTION

Risk of patient or user injuries which may be severe
Risk of damaging the product or environmental damage

PLEASE NOTE / IMPORTANT

Useful additional information and recommendations

PROHIBITED


AUTHORIZED

1.2. Description of the DolorClast® Focused Shock Waves


The DolorClast® Focused Shock Waves system (console, handpiece and gel pads) is designed to provide focused extracorporeal shock wave therapy (fESWT). A high-voltage pulse, controlled by the console, is used to excite piezoceramic elements located inside the handpiece arranged on a concave surface, causing them to expand briefly and simultaneously by a few micrometers, and generate a shock wave. Please see chapter "3. Components", page 8 Figures 1, 2 & 3 for device and part views.

1.3. Operating principle of the handpiece

The shape of the emission source (the handpiece) focuses the shock waves (S: 5 MPa penetration up to 65 mm) generated by the piezoceramic elements incorporated into the handpiece. Please refer to the following handpiece diagrams, which show the targeted midpoint focus at a depth of 40 mm (F).

The energy is delivered into the body in a focused manner, and the treatment depth varies in line with the type of gel pad used. The gel pads are available in the following three sizes, whereby E is the focus penetration depth and F is the constant midpoint focus. So, the thicker the gel pad, the shallower the penetration depth (E).

Gel pads	ID	E		Optimal depth of treatment (mm)**
		Midpoint focus penetration depth*	F (mm)	
Shallow	S	10	40	0–20
Medium	M	25	40	15–35
Deep	D	40	40	30–50

* Distance between the midfocus point and the edge of the gel pad.

** Depends on the target tissue depth and individual anatomical variants. Highest energy is localized at the midpoint focus, noting that at 5 MPa, the penetration depth achieved is up to 65 mm. See chapter "1.13. Indications", page 7.

1.4. Intended use

The DolorClast® Focused Shock Waves system is intended to be used for the treatment in the various fields indicated in chapter "1.5. Indications for use", provided that the affected parts of the body (tendons, muscles, bones, ligaments, and skin) can be safely reached and treated by focused extracorporeal shock waves.

1. PLEASE READ THIS FIRST (CONT.)

1.5. Indications for use

Field of Orthopedics

Acute or Chronic Tendinopathies:

- Tendinitis of the shoulder / Subacromial pain syndrome or impingement syndrome
- Achilles' tendinopathy or Achillodynia
- Lateral epicondylitis, i.e., tennis elbow, and radial/ulnar humeral/medial epicondylitis, i.e., golfer's elbow.

Musculoskeletal pathologies

- Greater trochanteric pain syndrome
- Chronic or acute muscular pain with no pathological correlate of injuries such as strains, bruises, sprains, muscle fiber torn, i.e.,:
 - Trigger points and myofascial pain syndrome
 - Stimulation of acupuncture points
 - Chronic lumbar and cervical pain with no pathological correlate of injuries and diseases such as vertebral fracture, herniated vertebral disk, spondylodiscitis, spondylolisthesis, ankylosing spondylitis, primary and secondary bone tumours, bacterial and/or viral inflammations.
- Plantar Fasciopathy
- Shoulder adhesive capsulitis or frozen shoulder

Bone pathologies

- Nonunion or delayed union fractures, stress fractures,
- Medial Tibial stress syndrome or tibial edge syndrome or tibial stress syndrome or tibial tendonitis
- Knee osteoarthritis

Field of Dermatology

- Chronic soft tissue wounds, including burn wounds, and diabetic foot ulcers

Field of Neurology

- Spasticity

Field of Urology

- Erectile Dysfunction
- Chronic pelvic pain syndrome

1.6. Intended clinical benefits

The intended clinical benefits of the DolorClast® Focused Shock Waves, is to be safely and effectively used for the treatment of the indications for use described in chapter "1.5. Indications for use".

1.7. Intended user and environment

The DolorClast® Focused Shock Waves system must only be used by qualified healthcare professionals (such as doctors and physiotherapists) in hospitals, clinics, medical universities, and private practices.

This product must be used in full compliance with the national regulations applicable in the country of use, and its use must also strictly respect accident prevention measures and these Instructions for Use.

No specific training is required to use this medical device.

The practitioner is responsible for performing the clinical treatments, and is liable for any risks or dangers that may arise due to lack of skill and/or medical expertise.

1.8. Intended patient population

The DolorClast® Focused Shock Waves system is intended to be used in the following target patient populations:

- Adult population (including breastfeeding women) suffering from the conditions set out in the "Indications" section, and who do not present any contraindications.

1.9. Contraindications

 The use of the product is contraindicated in the following conditions:

- Treatment over tissues filled with gas or air (namely the lungs or intestines),
- Treatment of patients under the age of 18,
- In the case of spasticity, treatment over growth plates is contraindicated,
- Pre-ruptured tendons, and fascias
- Pregnancy,
- Treatment over electronic medical implants such as pacemakers,
- Treatment of patients with blood-clotting disorders (including local thrombosis),
- Treatment of patients treated with oral anticoagulants and hemodiluting medication,
- Treatment of tissue presenting local tumors or local bacterial and/or viral infections,
- Treatment of patients under cortisone treatment (within the 6-week period following the administering of the last local cortisone injection).

1.10. Side-effects

The following potential side effects can occur from the use of DolorClast® Focused Shock Waves system. The side effects subside relatively quickly (after 1-2 weeks):

- Blood vessel injury due to misdirection of ESWT energy than can lead to local ecchymosis and bruising (mild superficial or subdermal hematomas), petechiae,
- Swelling and local oedemas,
- Pain and/or fatigue and/or discomfort during or after the shock wave treatment, and/or soreness or tenderness in the treated area (aggravation of pain scores/ frequency),
- Nerve injury due to misdirection of ESWT energy that can lead to mild neurological symptoms (numbness, tingling, throbbing, burning sensation),

1. PLEASE READ THIS FIRST (CONT.)

- Neurovegetative symptoms (nausea, dizziness, sweating, headache, tremors, fainting, dysesthesia),
- Local skin lesions (excoriation, blistering, redness, irritation),
- Dysuria when treating erectile dysfunction.

1.11. Potential complications

The following potential complications can occur from the use of DolorClast® Focused Shock Waves system only in the case of failure to respect the contraindications:

- Potential damage to the tissues in the lungs and intestines,
- Potential rupture in pre-ruptured tendons and fascias,
- Potential damage to immature joint cartilage epiphyses,
- Possible dissemination of malignant tumors,
- Skin lesions may develop if patients are pre-treated with cortisone for a longer period of time,
- In case of local anesthesia (not recommended): anesthesia complications, including allergic reactions to local or regional anesthetic agents,
- Propagation of existing infection (including fever reaction): Bacterial infection, Application site infection, Cellulitis, wound infection, Osteomyelitis.
- Transient heart rhythm disorder

1.12. Warnings and precautions

These Instructions for Use aim to ensure the correct installation and use of the DolorClast® Focused Shock Waves. Always keep these instructions close to hand.

The application settings outlined in this manual (handpiece, gel pad, number of sessions, number of impulses per session, frequency, etc.) are non-binding and are purely provided for recommendation purposes. They are based on reports provided by experts in Focused Extracorporeal Shock Wave Therapy (FESWT), published in the scientific literature, and presented at international conferences.

The patient is at the center of the treatment process, and their feedback is essential in locating the painful area and ensuring that the intensity is set just below the patient's pain threshold. We recommend starting treatment with low intensity settings.

Improper use of the DolorClast® Focused Shock Waves system may result in:

- Potential damage to the tissues in the lungs and intestines
- Potential rupture in pre-ruptured tendons and fascias
- Potential damage to immature epiphyses
- Possible dissemination of malignant tumors

- Potential damage of joint cartilage
- Skin lesions may develop if patients are pre-treated with cortisone for a longer period of time.

When treating patients with infectious (contagious) disease, the following measures must be taken:

- Before the treatment, cover any open wounds so as to ensure that no body fluids can leak out.
- Disinfect the device after the treatment. Follow the instructions provided by the manufacturer of the disinfectant product used.

Avoid applying shock waves to organs which contain trapped air (e.g. lungs, intestines). Double-check the entry window and position of the handpiece, and ensure the patient is in a suitable position. In addition, shock waves can cause adverse cardiac reactions.

Ultrasonic oscillations may negatively affect the function of cardiac pacemakers and ICDs (implanted cardio-defibrillators). Therefore, when using this product, we recommend that you do not treat patients over their cardiac pacemaker or ICD.

The efficacy of ESWT may be reduced in patients with hypogonadism.

The effects of ESWT on subsequent graft success in chronic wound patients are unknown.

Use for the intended purpose only. Before using the product, make sure that you have fully read and understood the operating instructions. This also applies to any equipment used alongside this product. Failure to strictly follow the operating instructions may result in serious injury or harm to the patient or user, or may result in the product being damaged, possibly beyond repair.

Risk of explosion: Do not use this product in the presence of inflammable anesthetics or gases.

This product must only be used by a physician or under their instruction.

Always inspect the product, checking for damage, before starting treatment. Damaged accessories or a damaged unit must not be used, and must be replaced. Only original EMS spare parts and accessories must be used.

Repairs should only be carried out by an EMS approved repair center.

1. PLEASE READ THIS FIRST (CONT.)

1.13. Indications

- The DolorClast® Focused Shock Waves system has been designed for use in various fields, such as the indications listed in the following table.
- The choice of the gel pad is under the user's responsibility.
- The choice of the gel pad depends on the target tissue depth and individual anatomical variants. The following table provides suggestions per indications.

⚠ For further information, refer to chapter "1.3. Operating principle of the handpiece", page 4 on how to estimate the optimal depth of treatment.

		Suggested gel pads		
		S	M	D
				
		Optimal depth of treatment (mm)*		
		0–20	15–35	30–50
INDICATIONS LIST				
Orthopedics				
Acute or Chronic Tendinopathy				
Shoulder tendinitis / Subacromial pain syndrome or shoulder impingement syndrome		✓	✓	
Achilles tendinopathy or Achillodynia		✓	✓	
Lateral epicondylitis, or tennis elbow, and radial/ulnar humeral/medial epicondylitis, or golfer's elbow.		✓		
Musculoskeletal pathologies				
Greater trochanteric pain syndrome			✓	✓
Chronic or acute muscular pain with no pathological correlation with injuries such as strains, bruises, sprains, muscle fiber tears, etc., or: Trigger points and myofascial pain syndrome		✓	✓	
Stimulation of acupuncture points		✓		
Chronic lumbar and cervical pain with no pathological correlation with injuries or diseases, such as vertebral fracture, herniated vertebral disk, spondylodiscitis, spondylolisthesis, ankylosing spondylitis, primary and secondary bone tumors, bacterial and/or viral inflammation.		✓		
Knee osteoarthritis				✓
Plantar Fasciopathy		✓	✓	
Bone pathologies				
Nonunion or delayed union fractures, stress fractures		✓	✓	✓
Medial tibial stress syndrome or tibial edge syndrome or tibial stress syndrome or tibial tendinitis		✓		
Exostoses, i.e., bone spur, Haglund's deformity		✓	✓	✓
Adhesive capsulitis or frozen shoulder			✓	✓
Dermatology				
Chronic soft tissue wounds, including burn wounds, and diabetic foot ulcers		✓	✓	
Neurology				
Spasticity		✓	✓	
Urology				
Erectile dysfunction		✓	✓	
Chronic pelvic pain syndrome		✓	✓	✓

2. WARNING

 EMS and the distributor of this product may not be held liable for any direct or consequential injury or damage resulting from its improper use, arising namely due to failure to follow the operating instructions, or incorrect or insufficient preparation and maintenance.

 Before using this product (including the console, the handpiece, and accessories), please take the time to carefully read, understand, and implement the recommendations set out in these Instructions for Use. Failure to observe the operating instructions may result in serious injury to the patient or user and/or in damage to the product. This product may only be used for its intended use by qualified professionals and for the applications described in this manual. If the product is used in combination with other instruments or devices, please refer to their respective instructions for use.

 The product must always be installed in such a way that permits the functioning of its ventilation. Never obstruct the back of the device (with a towel, etc.). It should be positioned at a minimum distance of 10 cm from walls.

 The product is disconnected from the mains power supply when the power cable connector is disconnected from the console. The power cable connector must be accessible at all times.


 Do not touch the patient and the Data Interface at the same time during the treatment.

 The device has been tested and certified as compliant with the standard IEC 60601-1 regarding electromagnetic compatibility restrictions. However, it is nevertheless possible that the device may cause weak electromagnetic interference with other devices or may be subject to minor interference from other equipment.

3. COMPONENTS

 The components delivered with your product will vary according to your configuration.


3.1. Figure 1: Unpacking

No.	NAME	QTY
1	"Instructions for Use" (IFU) manual	1
2	DolorClast® Compact Cart – Optional (screws included)	1
3	DolorClast® Focused Shock Waves console	1
4	Handpiece holder kit	1
5	DOLORCLAST® FOCUS handpiece	1
6	Gel pad kit: 3 gel pads + 1 storage box	1
7	Gel bottle for applying gel onto the gel pads	1
8	Power cable	1

3. COMPONENTS (CONT.)

3.2. Figure 2: Handpiece and gel pad kit


No.	NAME	QTY
1	Lock ring (available as a spare part)	1
2	Non-detachable handpiece cup: Receptacle for the gel pad	1
3	Minus [-] button: Decreases the intensity	1
4	Plus [+] button: Increases the intensity	1
5	ON/OFF button: Start, stop or pause the treatment	1

No.	NAME	QTY
6	Gel pad kit cover	1
7	Gel pad storage box	1
8	Deep* gel pad: D	1
9	Medium* gel pad: M	1
10	Shallow* gel pad: S	1

*See chapter "1.3. Operating principle of the handpiece", page 4.

3.3. Figure 3: Console

No.	NAME	QTY
1	Mains power switch	1
2	Touchscreen: Human Machine Interface. See chapter "3.4. Touchscreen interface introduction – Main HMI screens", page 10	1
3	Handpiece connector (applied part)	1
4	Red dot connector indicator for handpiece connection	1
5	Handpiece holder	1
6	Power supply socket (contains the fuse compartment)	1
7	Equipotential connector	1
8	Data interface (removable cover) – Reserved for AFSS (After-Sales Service)	1

No.	NAME	QTY
9	Product ID label: Serial number	1
10	Label: Information on the connector	1
11	Screw hole for securing the console to the DolorClast® Compact Cart	1

Please note that the data interface is only used for data transfer and cannot be used as an external power supply for other devices.

The data interface must only be connected to products with safe voltages and configurations that comply with Section 16 of IEC-60601-1, and must meet at least two MOOP with the mains power elements.

3. COMPONENTS (CONT.)

3.4. Touchscreen interface introduction – Main HMI screens

 This HMI is a touchscreen. The active menus or modes will appear highlighted on the display of the touchscreen.

Please note:

- There are 3 landing screens from where the user can select treatments. These treatments are to be selected in line with your best practices, and your knowledge of the patient and their condition.
- The main product functions that can be accessed via the touchscreen:
 - Registering your EMS product with a QR code.
 - Setting up the background, language, date and time, landing screen when switching on your product (Treatment, Indication, Custom).
 - Changing and saving certain preset protocols with the Custom function, so you can adapt the treatment while decreasing or increasing the energy, frequency, or impulses, or while using the +/- buttons on the handpiece.

+/- buttons on the handpiece.


- Using the special modes **Ramp up** or **Analgesic**, which help practitioners to improve sensitive patients' tolerance of the treatment.
- Selecting the function allowing **Ramp up** settings, in order to set your preferred scales for Intensity steps and Impulses
- Activating or deactivating the following functions enables:
 - **VAS assessment:** To display the scales relating to the patient threshold pain before and after treatment.
 - **Tutorials:** To view medical information, accessible from the Indication screen.
- Obtaining a summary of your treatment data, including the patient VAS assessment results (scale from 1 to 10).

3.4.1. The Ready-to-Work screens: Treatment, Indication, Custom


Please note that:

- The default screen displayed when you power the device on is the one you have selected as your landing screen. The default screen is the **Treatment** screen (as illustrated below). The user may choose a different landing screen between Indication screen and Custom screen by clicking on the icon and selecting the path  / **Console settings** / **Landing screen**.
- The blue buttons show your current selection.


Treatment screen

Indication screen



Custom screen

3. COMPONENTS (CONT.)

3.4.1. Interface introduction – Main HMI screens (cont.)

Treatment screen	
SELECTION	ACTION
	Press the icon to switch between INTENSITY and ENERGY (unit = mJ/mm ²).
	Pressing the – or + buttons will respectively decrease or increase the following main parameters: INTENSITY which can be switched to ENERGY To adjust the energy of the shock waves transmitted into the patient's tissue, increase or decrease the Intensity either via the Interface, using the + or – buttons, or by using the +/- buttons on the handpiece.
	FREQUENCY To adjust the number of impulses per second, increase or decrease the frequency. The frequency can be set to between 1 and 8 Hz, in increments of 1 Hz.
	PULSES To adjust the total number of impulses transmitted to the patient, increase or decrease the number of impulses. The number of impulses can be set to between 100 and 4,000, in increments of 100. You can also refer to chapter "6.3.1. Manual settings: Intensity / Energy, Frequency, Impulses", page 19.
	Displays the SETTINGS screen. Please refer to chapter "3.4.2. SETTINGS Screens: Console - Treatment", page 12.
Preset	Displays the Indication / Custom screens.
Total energy X/YYY mJ/mm ²	Cumulated energy density in mJ/mm ² delivered since the beginning of the treatment.
Total duration XX min YY s	Total cumulated duration in minutes and seconds that have elapsed since the beginning of the treatment.
	<u>Special mode*</u> : 1. If the Analgesic mode is enabled, the minimum intensity level is reduced to 0.1, and can be increased up to 1, by steps of 0.1. From intensity level 1, the intensity is increased by steps of 1 up to the maximum intensity of level 20. See also "Special modes: Analgesic and Ramp up", page 16. 2. If the Analgesic mode is disabled, the intensity minimum value is 1. The intensity is increased by steps of 1 up to the maximum intensity of level 20
	<u>Special mode*</u> : This function automatically adds or subtracts Intensity steps and Pulses for the selected treatment, with a percentage of energy below or above the minimum setting. If Ramp up number 1 is selected, the intensity will automatically increase by 1 unit every 150 impulses (no further manual action is needed). This preset increase cannot be changed by the user. If Ramp up number 2 is selected, then your current customized Ramp up 2 version will be enabled. You can define up to eight increment settings for the intensity steps and impulses. Please refer to "Special modes: Analgesic and Ramp up", page 16. You can add or remove certain increasing or decreasing scales from the SETTINGS screen / Treatment settings / Configure Ramp up 2 settings .

- * The special modes can only be accessed via the **Treatment** screen or with the  icon via any other screen.
* Only one special mode can be selected at a time.

Indication and Custom screens

SELECTION	ACTION
	Takes you back to the Treatment screen.
Indication**	Displays the preset protocols, which list different indications with recommended settings and photos. If the Tutorials function is activated via the SETTINGS screen / Treatment settings , that will enable these Tutorials to be displayed. See chapter "5.2.4. VAS assessment and Tutorials", page 17.
Custom**	Displays your pre-defined Protocols. See chapter "6.3.4. Creating or Saving your custom protocol", page 21.

3. COMPONENTS (CONT.)

3.4.2. SETTINGS Screens: Console - Treatment

The blue buttons show your current selection.

SETTINGS – Main buttons	
SELECTION	ACTION
	Takes you back to the Treatment screen
	Takes you back to the previous screen

3. COMPONENTS (CONT.)

3.4.2. SETTINGS Screens: Console - Treatment (cont.)

Console settings	
SELECTION	ACTION
1. Display	1. Allows you to change the background color: Dark or Light . The blue buttons show your current selection.
2. Brightness	2. Allows you to change the background brightness. The blue buttons show your current selection.
3. Standby	3. Enables the touchscreen standby mode. chapter "5.2.2. General modes", page 16.
4. Language	4. Allows you to select your language. The button on the right displays your current selection.
5. Date/time	5. Allows you to select the current date. The button on the right displays your current date/time.
6. Landing screen	6. Allows you to choose between Indication and Custom for the first landing screen, instead of the Treatment screen, which will load when you switch on the console. The button on the right displays your current selection.
7. System logs	7. Allows you to save device log information from the product onto a USB stick, for After-Sales Service purposes.
8. Commissioning/ Next maintenance	8. Set the commissioning date and the date of the next maintenance due.
9. After-Sales	9. This option is reserved for After-Sales Service team members only.
10. Register	10. Enables the console to be registered via a QR code.

Treatment settings	
SELECTION	ACTION
1. VAS assessment	1. You can activate or deactivate the VAS assessment function at any time: If it is activated, you will be able to display the patient pain assessment. chapter "5.2.4. VAS assessment and Tutorials", page 17.
2. Tutorials	2. If activated, this will mean that you can access the Tutorials screen via the Indication screen. chapter "5.2.4. VAS assessment and Tutorials", page 17.
3. Ramp up 2 settings	3. You can select Ramp up mode numbers 1 or 2 from your setting screens. Ramp up 1 is the default mode, and Ramp up 2 needs to be custom configured before use. See "Special modes: Analgesic and Ramp up", page 16.

3.4.3. Factory default settings

PARAMETER	DEFAULT VALUE	RANGE	STEP OF
Intensity	5	0.1–20	0.1 or 1
Energy	0.064 mJ/mm ²	0.008–0.403 mJ/mm ²	0.001–0.048 mJ/mm ²
Number of impulses	2,000	100–4000	100
Frequency	6 Hz	1–8 Hz	1 Hz
Analgesic mode	OFF	OFF, ON	N/A
Ramp up mode**	OFF	OFF, ON*	N/A
Tutorials**	On	Off, On	N/A
VAS assessment**	On	Off, On	N/A

* Intensity = 1 every contact when selecting Ramp up mode

** To be set with Treatment settings

4. INSTALLATION

4.1. Installing the console

⚠ Install the console on a flat, stable surface or use the cart (optional) designed specifically for the console.


⚠ The cart has been designed exclusively for the use of compatible EMS products.

Please only install the DolorClast® Focused Shock Waves device and related accessories on this cart, as indicated in Chapters 4.1.1 and 4.1.2, and no other equipment.


4.1.1. Installing the console on the cart

All screws must be pre-assembled, without tightening them at first. They should only be tightened once they are all in the correct positions.


1. Make sure that the brakes on the cart are engaged before installing and using the device.

2. Place the console on the cart (if applicable). Make sure that the feet are correctly aligned with the corresponding notches on the upper surface of the cart.


3. Secure the console to the cart using the dedicated thumbscrew.


4.1.2. Installing the handpiece holder onto the device

1. Use the screws and tool provided to secure the handpiece holder onto the console.

2. Make sure that the handpiece holder is in the correct position. Please refer to Section 5 of chapter "3.3. Figure 3: Console", page 9.


4. INSTALLATION (CONT.)

4.2. Connecting the console to the equipotential conductor

-  When applicable, and according to your internal protocols, connect the equipotential connector to the rear of the console, with the busbar.
-  The equipotential cable is not supplied with the console.


The equipotential connector is used to prevent differences of potential between items of medical equipment, which could generate errors in measurements or treatments, if the medical device generates interference.

If interference is observed, the unit must be connected to other devices via this plug.


4.3. Connecting the handpiece to the console

1. Connect the handpiece to the console. Take care to ensure that the connector aligns with the red dot.


-  Plugging and unplugging the handpiece must be done exclusively holding the metallic connector. Do not attempt to plug or unplug the handpiece by pushing or pulling the cord.
- 2. Position the handpiece in the handpiece holder. The handpiece cable can be wrapped around the handpiece holder.

4.4. Connecting the power cable

To prevent the risk of electric shock, this product must only be connected to a mains power supply which is fitted with protective earthing.

Connect the power cable to the power socket located at the rear of the console. See Section 6 of chapter “3.3. Figure 3: Console”, page 9.


5. GETTING STARTED

5.1. Starting up the console

- ⚠ The product's mains power switch must be accessible at all times, so that the device can be easily disconnected, if needed.
- ⚠ Make sure that the handpiece is connected.

Press the mains power switch, which is located on the left-side panel.

5.2. Interface

- ⚠ In case of any issues during start-up, please refer to the Troubleshooting guide chapter "9. Troubleshooting - Warranty - Support", page 24. Examples: Error message, fuse change required.
- ⚠ If the software, translations or indications have to be updated, please contact your after-sales service. Please refer to the troubleshooting guide.
- ⚠ Once the console and handpiece are connected, your **Landing screen** is displayed, giving you easy access to the **Indication**, **Custom**, and **Treatment** screens. Please refer to chapter "3.4.1. The Ready-to-Work screens: Treatment, Indication, Custom", page 10.

5.2.1. Configure your background

From the **Treatment** screen / Icon  / **SETTINGS** / **Console settings**.

Display and Brightness.

The default initial background screen can be brightened by sliding the icons in these modes.

5.2.2. General modes

Standby

If enabled: If the device isn't used for 30 minutes, the touchscreen will automatically go into standby mode. Standby mode is disabled when you touch the screen or press one of the buttons on the handpiece.

Language, Date/Time

Press the corresponding icons to update this information, if needed.

Landing screen

Once selected, choose the first screen to be displayed when you start up the console.

5.2.3. Special modes: Analgesic and Ramp up

- 👉 Only one special mode can be selected at a time. These special modes can be accessed during treatment and from the **Treatment** screen, and via the home icon: , you can select one of the two special modes: **Ramp up** or **Analgesic**. See chapter "3.4.1. The Ready-to-Work screens: Treatment, Indication, Custom", page 10 and the table descriptions on the following page.

Analgesic mode

The **Analgesic** mode helps practitioners to improve sensitive patients' tolerance of the treatment.

For further information on this special mode, please refer to the Treatment screen table in chapter "3.4.2. SETTINGS Screens: Console - Treatment", page 12.

Ramp up mode

The **Ramp up** modes are designed to help you treat a patient at maximum intensity, without having to manually adjust the intensity. Shock wave treatments are dose-dependent. The best treatment outcomes will be achieved at the maximum possible intensity tolerated by the patient.


If the treatment is paused, the **Ramp up** mode is disabled. Press the ON/OFF button to continue the treatment, using the last intensity reached.

5. GETTING STARTED (CONT.)


Ramp up 1 is the default mode, and increases the intensity by one unit after every 150 impulses.

Ramp up 2 can be set:

1. From the **Treatment** screen / Icon  to **SETTINGS / Treatment settings** / select **Configure** the **Ramp up 2 settings**.

2. The intensity and Impulse scales can be increased or decreased using the - and + buttons.
- Once you **Remove** or **Add**, the **Preview** is updated. The **Ramp up 2** settings allow you to add up to 8 steps of Intensity / Impulses.

3. Save to go back to your previous screen.



5.2.4. VAS assessment and Tutorials

You can activate or deactivate the **VAS assessment** and **Tutorials** functions at any time:

From **Treatment** screen / Icon  / **SETTINGS / Treatment settings** / select **On** or **Off** for one or both of the following functions:

- **VAS assessment:** If activated, a VAS assessment screen will be displayed before and after the treatment. This screen will allow the patient to evaluate their pain level directly on the touchscreen. See chapter “6.4.3. Assessing patient pain”, page 22.

- **Tutorials:** If activated, this will enable the **(Tutorials)** screen to be displayed just after clicking on the **Ready** button, which can be accessed from the **Indication** screen once you have selected a protocol, for example:

6. TREATMENT

⚠ Always store the handpiece in its holder, which is mounted on the side of the console: before and after treatment, and when pausing.

👉 If you have enabled the **Stand by** function, after 30 minutes without use, the touchscreen will automatically go into standby mode. Standby mode is disabled when you touch the screen or press one of the buttons on the handpiece. Please refer to “3.4.2. SETTINGS Screens: Console - Treatment”.

6.1. Prerequisites

⚠ The practitioner must be thoroughly familiar with the product and the first chapters of these Instructions for Use.

⚠ Take care when using excessively high energy shock waves. High energy shock waves may cause tissue damage. Do not apply a higher number of shock waves during one treatment session than the number specified.

⚠ The product must be checked, cleaned and tested according to chapter “5. Getting started”, page 16 and chapter “7. Reprocessing and Maintenance”, page 23.

6.2. Preparing the patient and the handpiece

⚠ Use only EMS gel pads for patient treatment. Refer to chapter “1.13. Indications”, page 7 to identify the appropriate gel pad for the treatment you'll be administering. Please refer to chapter “6.3.3. Factory settings: Indications and preset protocols”, page 20. The preset protocols are selected according to the location of the treated area, and depending on slight individual skin variations.

⚠ Once you have applied the contact gel, there should be no air bubbles between the handpiece and the gel pad, or between the gel pad and the patient's skin, in order to ensure the effective and correct transmission of the shock waves, otherwise you risk losing treatment efficiency.


• For optimal results, use the EMS SWISS DolorClast® contact gel. Thanks to its special viscosity, the EMS gel transmits shock waves without any leaks or diffusion, ensuring perfect transmission throughout treatment. Do not change the gel pad while the handpiece is operating.

👉 The gel pads are delivered ready-to-use. For further information on gel pad maintenance, please refer to the “Maintenance” chapter.

👉 The number of sessions is to be determined at the practitioner's discretion.

1. Get the patient ready in the treatment room.
2. Take the patient's medical history in order to determine any medications the patient may be taking and any other treatments which they may be under.
3. Establish a diagnosis of the pathology, or check the diagnosis received.
4. Rule out any contraindications.
5. Explain to the patient how shock wave treatment works: The kinetic energy created by pressure waves is transferred to the handpiece, and then converted into a shock wave, which is delivered straight to the injured tissue. It accelerates tissue repair and relieves pain.
6. Select the appropriate gel pad and connect it to the handpiece, as follows,
7. Position the patient on the treatment bed.

Figure 4 - Preparing the handpiece and gel pad

- A. Using the padlock markings as a guide, unscrew the lock ring.
- B. Apply a sufficient amount of contact gel to the center of the cupped surface of the handpiece.
- C. Avoiding the development of air bubbles, position the rounded part of the gel pad against the cupped surface of the handpiece, and use a twisting motion to spread the gel evenly.
- D. Align the markings to screw the lock ring back into place.

👉 If the lock ring is incorrectly positioned, this will cause the handpiece to malfunction. Make sure that the handpiece's lock ring is fully tightened.

6. TREATMENT (CONT.)

8. Use palpation to locate the area to be treated.
Use palpation to find the most painful area.

-  Use biofeedback from the patient to locate the affected tissue.

9. Once you have located it, mark the skin to clearly identify the treatment area.

10. Once the patient is in position and you have identified the area to be treated, apply a small amount of gel to the treatment area.

11. Use the gel pad to spread the gel onto the treatment area.


6.3. Preset protocols and treatment settings

Generally, the treatment settings include the settings for Intensity/Energy, Frequency, Impulse.

- The preset protocols are available in the Indication function to help you target the treatment and increase efficiency.
- At any time, you can adapt and save your own preset protocols via **Custom** and **Indication**.
- The displayed lists will include the Indications factory settings and your saved treatments.

6.3.1. Manual settings: Intensity / Energy, Frequency, Impulses


The **INTENSITY** or **ENERGY**, the **FREQUENCY** and the number of **IMPULSES** can be adjusted by pressing the - or + icons.


Please refer to the **Treatment** screen table in chapter “3.4.2. SETTINGS Screens: Console - Treatment”, page 12.

6.3.2. Indication screen

- At any time, the **Tutorials** can be deactivated via the **SETTINGS** screen. See chapter “5.2.4. VAS assessment and Tutorials”, page 17.
 - Deactivating the **Tutorials** will deactivate their display before starting your treatment.
1. Via the **Indication** screen, select your concerned preset protocol will lead to the factory settings and **Medical information**.

2. Selecting **Medical information** will display the general medical information.

3. If activated, selecting the **Ready** button will add the **Tutorials** before starting the treatment, otherwise the preset protocols will be displayed directly.
4. From the next step, you can still edit the default settings. Click on **Start** when you are ready to begin the treatment.

6. TREATMENT (CONT.)

- These are examples which depend on the location of the treatment area, and small individual skin variations.
- The number of sessions is to be determined at the practitioner's discretion.
- The choice of gel pad depends on the target tissue depth and individual anatomical variants. Suggestions are provided to guide you in selecting the type of gel pad to use.

6.3.3. Factory settings: Indications and preset protocols

The preset protocols are determined in regard to the average parameters, as described in clinical references.

INDICATIONS - PRESET PROTOCOLS	Factory settings		
	Intensity / energy	Frequency (Hz)	No. of impulses
Orthopedics			
Acute or Chronic Tendinopathy			
Shoulder tendinitis / Subacromial pain syndrome or shoulder impingement syndrome	12	5	1,500
Achilles tendinopathy or Achillodynia	5	5	1,500
Lateral epicondylitis, or tennis elbow, and radial/ulnar humeral/medial epicondylitis, or golfer's elbow.	5	5	2000
Musculoskeletal pathologies			
Greater trochanteric pain syndrome	5	5	1,500
Chronic or acute muscular pain with no pathological correlation with injuries such as strains, bruises, sprains, muscle fiber tears, i.e.			
Trigger points and myofascial pain syndrome	5	5	2,000
Stimulation of acupuncture points	5	5	2,000
Chronic lumbar and cervical pain with no pathological correlation with injuries or diseases, such as vertebral fracture, herniated vertebral disk, spondylodiscitis, spondylolisthesis, ankylosing spondylitis, primary and secondary bone tumors, bacterial and/or viral inflammation.	5	5	2,000
Knee osteoarthritis	15	5	2,000
Plantar Fasciopathy	5	5	1,500
Bone pathologies			
Nonunion or delayed union fractures, stress fractures	17	3	2,000
Medial tibial stress syndrome or tibial edge syndrome or tibial stress syndrome or tibial tendinitis	5	5	1,500
Exostoses, i.e., bone spur, Haglund's deformity	5	3	2,000
Adhesive capsulitis or frozen shoulder	8	5	1,200
Dermatology			
Chronic soft tissue wounds, including burn wounds, and diabetic foot ulcers	5	5	2,500
Neurology			
Spasticity	5	5	2,000
Urology			
Erectile Dysfunction	5	5	2,000
Chronic pelvic pain syndrome	15	3	3,000


6. TREATMENT (CONT.)

6.3.4. Creating or Saving your custom protocol


 You have the option of adapting your settings, as described in the following chapters.

Once you have adapted your personal preset protocols, by adjusting the Intensity, Frequency, and Number of Impulses, you will be able to save or delete your custom settings, as needed.


1. If your landing screen is **Treatment**, select the **Preset** icon
2. Select the **Custom** screen


3. Either:
 - **Create a new custom**, which takes you straight to the Treatment creation screen. In the top-left of the screen, you'll see **Treatment creation**.
 - Alternatively, you can select your custom treatment which you have saved previously. In the top-left of the screen, you'll see the name of your selected treatment type.
4. Adjust your settings if needed. Please refer to the **Treatment** screen table in chapter "3.4.2. SETTINGS Screens: Console - Treatment", page 12 for details on adjusting the **INTENSITY** / **ENERGY**, **FREQUENCY**, **IMPULSE**

5. **Save** your custom treatment. Please note that you can **Delete** your previously-saved custom treatment.
6. Define and confirm the **Body location** or press the **Skip** button.

7. Check or correct the proposed name (up to 22 characters).

8. Press save, and you'll be taken back to your list of custom treatments.

6.4. Performing the treatment

 Do not use the treatment directly on cracked or broken skin. If treatment is necessary, use a treatment interface, such as an adhesive waterproof polyurethane film (for example, the OPSITE™ wound dressing from Smith and Nephew or an equivalent) to protect the wounded area from direct contact with the gel pad.

 If your hands are getting tired, take a break before restarting treatment.

 Ensure the handpiece is cleaned and disinfected prior to treating a patient.

6.4.1. Handling the handpiece

 To maximize treatment efficiency, EMS recommends holding the handpiece in one of three main ways, always keeping the handpiece perpendicular to the treatment area.

- **A:** One-handed high position, enabling thumb activation of the buttons.
- **B:** One-handed low position, enabling positioning closer to the patient.
- **C:** Two-handed position, combining the two previous positions.

A

B

C

6. TREATMENT (CONT.)

6.4.2. Starting the treatment

1. Once the user has defined the settings on the Treatment screen, press **Start**.
2. If a preset protocol has been defined, press **Ready** to access the following screen and once on the Treatment screen press **Start**.
3. If your settings are ready (from **Treatment** or **Custom** or **Indication** screens), select **Start** and click on the handpiece button to start treatment.

 The application force exerted manually varies, and it depends on the indication. See chapter “1.13. Indications”, page 7.


 For patients with a very low pain threshold, treatment can be started with the **Analgesic** mode or a **Ramp up** mode.

 Please refer to the “Special Modes” explained in chapter “3.4. Touchscreen interface introduction – Main HMI screens”, page 10.

6.4.3. Assessing patient pain

 This step can be carried out:
- Only if the **VAS assessment** function is enabled.
For more information, see chapter “5.2.4. VAS assessment and Tutorials”, page 17.

1. Ask the patient to answer the question displayed on the screen: **How severe is your pain today?**
2. The cursor must be positioned to indicate the patient's pain score.

3. Click on **Confirm**, and the score will appear on the screen.

 The Score displayed will also be included in the Treatment summary at the end of the treatment.

4. Select **Start** on the following screen.


6.4.4. Pausing / Ending the treatment

 If standby mode is enabled, the touch screen will automatically go into standby mode after 30 minutes of non-use. Standby mode is disabled when you touch the screen or press one of the buttons on the handpiece.

 Displaying your patient pain score is only possible if the VAS assessment mode is enabled.

 You can press the ON/OFF button on the handpiece at any time to pause the treatment.

1. Once the number of impulses remaining reaches “0”, select **End treatment** and the following screen will be displayed:
Treatment summary.


2. Place the handpiece back into the handpiece holder.
3. Wipe any residual gel off the patient.
4. Switch off the device.
5. Disconnect the handpiece from the unit.
6. Proceed with the maintenance steps, as described in the following chapter.

7. REPROCESSING AND MAINTENANCE

Daily checks	Monthly checks
* Check the handpiece and the accessories for any damage, loose or missing parts, cleanliness, and completeness. * Check the handpiece connection cable for any possible damage.	* Any lettering, labeling or identification necessary for the safe intended use must be legible. * Any missing or illegible lettering, labeling or identification information which may lead to improper handling and reprocessing, must be restored. Please contact an EMS approved service center. * Check the contact surface of the handpiece for damage. * Check the power output of the handpiece.

7.1. Cleaning and disinfecting

- ⚠ Danger of electric shock! Make sure that no humidity gets into the handpiece. Do not immerse the handpiece or plug.
- ⚠ For reprocessing, disconnect the handpiece from the console. Make sure that no humidity gets into the handpiece. Do not use any cleaning products, abrasive agents or solvents on the handpiece.
- ⚠ Make sure that no humidity gets into the console. Do not use any cleaning products, abrasive agents or solvents on it.
- ⚠ EMS recommends cleaning and disinfecting the device using a disinfectant wipe with quaternary ammonium as the active ingredient (phenol-free and bleach-free), such as the CaviWipes™ from Metrex. DO NOT use the CaviWipes™ 2.0 from Metrex as they can damage our product.
- ⚠ Do not spray the handpiece with water. EMS recommends that you clean it with disinfectant wipes.

FREQUENCIES FOR REPROCESSING

Handpiece - Lock ring - Gel pads

Before and after each treatment.	1. Disconnect the handpiece from the console. 2. Unscrew the lock ring and remove the gel pad from the handpiece. 3. Clean the gel pad and the lock ring with warm water (Max. 30 °C), then clean the handpiece, gel pad and lock ring with a disinfectant wipe. 4. Reassemble the handpiece.
----------------------------------	--

Console

At the end of the day	1. Turn off the console. 2. Disconnect the power supply connector. 3. Clean and disinfect the surfaces using a disinfectant wipe.
-----------------------	---


- 👉 Always store the handpiece in its holder before and after use. Please refer to Step 2 of chapter “4.3. Connecting the handpiece to the console”, page 15.

7.2. Maintenance

- ⚠ To comply with the standard IEC-60601-1 relating to technical and safety tests, as well as maintenance intervals, the user must follow all applicable local laws and regulations in the country where the device is used.

If the fuses blow:

- ⚠ Replace exclusively with the type of fuse indicated on the rear of the unit. See Section 6 of “Figure 3: Console”, page 9.

1. Disconnect the power supply connector.
 2. Using a small screwdriver, remove the cover: Only use 250 V x 4 A fuses.
 3. Replace both fuses before putting the cover back on.
- 👉 This maintenance can be carried out by EMS or by an approved repair center. The fuses to be used are listed in chapter “11. Technical description”, page 27.
 - 👉 Should legal provisions in your country specify maintenance intervals, these must be observed. The console and handpiece may need to be returned for periodic servicing.
 - 👉 To avoid any incidents or damage caused by aging and wear, the product and its accessories should be serviced at appropriate intervals.
 - 👉 Depending on the frequency of use, but in any case, at least once a year, have an expert check the functional and operational safety of the equipment.

8. PRODUCT DISPOSAL

Product disposal: The product must not be discarded with domestic household waste.

Should you wish to definitively dispose of the product, please comply with the applicable regulations in your country.

Keep the original packaging until the product is to be disposed of permanently.

Waste Electrical and Electronic Equipment belonging to customers located in the European Union may be shipped to EMS for recycling in accordance with the WEEE regulations. The costs of recycling, exclusive of shipping fees, are covered by EMS.

9. TROUBLESHOOTING - WARRANTY - SUPPORT

⚠ CAUTION

- ⚠ Do not attempt to modify or repair the product yourself. Please contact an EMS approved service center.
- ⚠ The RTC battery cannot be replaced. Do not attempt to remove or replace the battery yourself.
- ⚠ Please make sure that your concerned troubleshooting case has been resolved. If some of the following issues persist, please contact the AfterSales Services. In such a case, you'll need the serial number, which can be found on a label positioned on the rear of the console. See Section 11 (Serial number) of chapter "3.3. Figure 3: Console", page 9. The serial number can also be revealed on the HMI, by following the path **SETTINGS / Console settings / Maintenance settings**.
- ⚠ If you encounter any issues with the product, please return all your previously disinfected devices (console, handpiece, gel pads) to the After-Sales Services.

Please refer to the following chapters regarding Technical Support and the Warranty.

- ⚠ Repair can be refused for products or accessories that are received in a contaminated condition.
- ⚠ The product must be tested before each treatment.
- ⚠ The console must be positioned horizontally.
- ⚠ Resetting the device will delete all custom and log data.
- 👉 As the log data is registered and in order to accelerate your After-Sales Service Support, your EMS service center may request this procedure:
Plug a USB2 stick, with 2-4 GB of capacity, format FAT32, into the rear of the console.
From the **Treatment** screen, select the **SETTINGS / Export** / Please follow the procedure displayed on the screen / **Export database**.

9.1. Troubleshooting table

DEVICE	CAUSE	POSSIBLE SOLUTION**
Console	The console is connected but there is no power at all	• Check your connections. • Check your general electrical power supply. • If the general electrical power supply is working normally, the internal console fuse in the device needs to be replaced. Please see chapter "7. Reprocessing and Maintenance", page 23.
Touchscreen	You no longer need your saved custom and your last saved Ramp up 2 settings	If you want the HMI to be reset to the default factory settings: From the screen Console settings / select Maintenance settings / select Reset . This will delete your Custom and Ramp up settings and the Analgesic mode and Tutorials mode, which will return to default settings.
Handpiece	1. ER1 code error is displayed on the HMI 2. Connected to the working console but it does not seem to be powered at all, or a warning screen is displayed on the HMI 3. Accidentally got wet	1. Check the handpiece connection. Please refer to the following error codes table. 2. Check the output power of the handpiece: press the handpiece power button while starting a treatment from the HMI. In case of displayed ER9 error code on your touchscreen, please refer to the following error codes table. 3. Stop this action as soon as possible. This error may cause serious damages to your device. Leave it to dry out, and refer to the following chapters and/or your After-Sales Service.

** As the log data is registered, and in order to accelerate the resolution of your issue with After-Sales Service Support, your EMS service center may request this procedure:
1- Plug a USB2 stick, 2-4 GB capacity, format FAT32, into the rear of the console.
2- From the **Treatment** screen, select the **SETTINGS / Maintenance settings / Export** / Please follow the procedure displayed on the screen / **Export database**.

9. TROUBLESHOOTING - WARRANTY - SUPPORT (CONT.)

9.2. Error codes

ERROR CODES	DEVICE	CAUSE	POSSIBLE SOLUTION***
ER2 ER3 ER4	Console	<i>High voltage fault</i> Possible failure** displayed message.	1. Restart the device. 2. Please contact your service center if the error persists.
ER5 ER6	Console	<i>Fan fault*</i> Possible warning* displayed message.	1. Click the "OK" button on the warning screen, and you can continue to use the device, as the console is equipped with 2 fans. 2. Restart the device. 3. Please contact your service center as soon as possible.
ER8 ER12	Console	<i>Console detection fault or internal error</i> Possible failure** displayed message.	1. Restart the device. 2. Please contact your service center if the error persists.
ER7	HMI	<i>HMI fault</i> Possible failure** displayed message.	1. Restart the device. 2. Please contact your service center if the error persists.
ER1	HP	<i>The handpiece is not detected</i> Possible failure** displayed message.	1. Connect your handpiece. 2. Restart the device if the error persists. 3. Please contact your service center if the error persists.
ER9	HP	<i>Handpiece fault</i> Possible warning* displayed message.	The handpiece may have reached its maximum recommended usage 1. Replace the handpiece. 2. Click the "OK" button on the warning screen, and you can keep on using the device with another handpiece, as long as you've connected it to the console.

* Warnings: Can be displayed on the touchscreen as long as the console starts, but you don't modify the device behavior. The error is registered in the log file.

** Failures: Can be displayed on the touchscreen as long as the console starts and as the energy delivery is still in progress. When the error occurs, the device might be shut down. The error is registered in the log file.

*** As the log data is registered, and in order to accelerate the resolution of your issue with After-Sales Service Support, your EMS service center may request this procedure:

1- Plug a USB2 stick, 2-4 GB capacity, format FAT32, into the rear of the console.

2- From the **Treatment** screen, select the **SETTINGS / Maintenance settings / Export** / Please follow the procedure displayed on the screen / **Export database**.

9.3. Warranty

The warranty is deemed null and void if the devices have been used with non-original EMS parts, accessories and handpieces. The warranty is void if the devices have been opened up/dismantled.

EMS and the distributor of this device may not be held liable for any direct or consequential injury or damage resulting from improper use of the device, namely caused by failure to observe the instructions for use, or by incorrect or insufficient preparation and maintenance.

9.4. EMS Technical Support

 Make sure that the device and its accessories have been used in accordance with the conditions specified by EMS.

 To accelerate the resolution of your issue and to improve the quality of our service, please provide the following information:

- Product reference number.
- Batch number/serial number, Service history of the device (e.g., previous issues or repairs).

Please contact your EMS approved service center for any product servicing or repairs. You must complete the appropriate EMS form in order to be issued with an Equipment Return Agreement (ERA) number.

EMS declines any responsibility for the safety of the device and declares the warranty invalid in the event of servicing or repairs being carried out by an unauthorized third parties, or in the event of non-original spare parts being used.

It is mandatory to return your product in its original packaging. By following these packaging guidelines, your product will be protected against damage in transport. To protect the staff at the EMS approved service center, and for reasons of safety in transport and shipment, all products and accessories returned to the factory for repair or servicing must be cleaned and disinfected in accordance with the instructions for use.

 Repair may be refused for products or accessories which are received in a contaminated condition.

When sending your product directly to the EMS authorized service center, please include the name of the distributor to simplify processing.

10. SYMBOLS / GLOSSARY

10.1. Symbols

SYMBOL	MEANING
	Manufacturer logo
	DolorClast® Focused Shock Waves
	Manufacturer
	Catalog number
	Serial number
	Date of manufacture
	Input power
	Risk of electric shock
	CE symbol refers to REGULATION (EU) 2017/745 and the notified body identification number
	Device requiring protective earthing
	Read the instructions for use
	Disposal of old electrical & electronic equipment (applicable in the European Union and other European countries with separate collection systems)
	Expiry date
	Maximum humidity of 85%
	Transport temperatures between +5 °C and +38 °C
	Medical Device

SYMBOL	MEANING
	Authorized representative in the European community
	CE marking
	EMS website
	Fragile
	Keep out of rain
	Vertical position
	Applied part, type B
	Equipotential plug
	Fuse
	Restriction of hazardous substances in electrical and electronic equipment
	Data transfer
	Unique Device Identification
	be dispensed only by authorized prescription of any physician or legally authorized health care provider under applicable state law
	Batch code
	Recyclable
	No stepping on surface

10.2. Glossary

ACRONYM	MEANING
VAS	Visual Analog Scale (used for the VAS assessment function)
FESWT	Focused Extracorporeal Shock Wave Therapy
HMI	Human Machine Interface
QTY	Quantity
AFSS	AFter-Sales Service
CONT.	Continued (for chapters headings)
IFU	Instructions For Use
CPU	Central Processing Unit
RTC	Real-Time Clock

11. TECHNICAL DESCRIPTION

MANUFACTURER		E.M.S. Electro Medical Systems S.A. Chemin de la Vuarpillière 31, 1260 Nyon, Switzerland	
PRODUCT NAME		DOLORCLAST® Focused Shock Waves	
MODEL		FT-249	
OPERATING CONDITIONS		+10 °C to + 35 °C Relative humidity: 35% to 65% Atmospheric pressure: 700 hPa to 1,060 hPa - Max. altitude: 3,000 m	
STORAGE AND TRANSPORT CONDITIONS		+5 °C to +38 °C Relative humidity: Maximum 85% Atmospheric pressure: 700 hPa to 1060 hPa	
REGULATION (EU) 2017/745 CLASSIFICATION		Console, Handpiece: Class IIa Contact Gel, Gel Pads, Compact Cart DolorClast®: Class I	
RANGES		Intensity: 0.1–20 Energy: 0.008–0.403 mJ/mm² Number of impulses: 0–4,000 Frequency: 1–8 Hz	
CONSOLE		HANDPIECE AND GEL PADS	
Classification EN 60601-1	Class I Applied part, Type B Mode: Continuous operation	Classification EN 60601-1	Class I Applied part, Type B
IP Classification IEC 60601-1	IP21	IP Classification IEC 60601-1	IP20
Supply voltage	100–127/220–240 VAC	Aperture diameter	100 mm
Rated power	150 VA	Aperture angle	94°
Mains frequency	50-60 Hz	Energy density (positive)	0.008–0.403 mJ/mm²
Noise level	70 dB(A)	Pressure range	6.4 to 80 MPa
Operating mode	Multiple impulse mode repetition rate: 1–8 Hz	Operating mode	Two +/- buttons: Decrease and Increase the Intensity One stop/pause button
Weight	Approx. 15 kg (approx. 33 lbs) Cart (optional): 18 kg (approx. 40 lbs)	Weight	Handpiece: Approx. 1.4 kg (approx. 3.1 lbs) Gel pad kit: Approx. 0.8 kg (approx. 1.8 lbs)
Dimensions in mm (H x W x D)	Console: 180 x 340 x 450 Cart (optional): 913 x 460 x 541	Dimensions in mm (H x W x D)	Handpiece: 156 x 116 x 116 Gel pad box: 380.4 x 140.4 x 80.2 Pads thickness H from focal base: 10-25-40
Service life	7 years	Service life	2 years or 5,000,000 pulses
Fuses	Two fuses: F4AL250V Interrupting rating: 63 A	Focal volume Intensity stages Overall energy E (5 mm) Optimal treatment depth 5 MPa penetration depth	142 mm³ 29 steps 23.58 mJ 5 to 40 mm 65 mm
		Focal size (min. to Max.): fx(-6dB) = fy(-6dB) fz(-6dB)	0.5 to 4 mm 3.5 to 17 mm

12. ELECTROMAGNETIC COMPATIBILITY

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

- The DolorClast® Focused Shock Waves system is intended for use in the electromagnetic environment specified below. The customer or the user of the DolorClast® Focused Shock Waves system should ensure that it is used in such an environment.
- These guidelines may not apply to all situations. Electromagnetic propagation is affected by absorption and reflection from surrounding structures, objects, and people. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radios, AM and FM radio broadcast, and TV broadcast cannot be theoretically predicted with accuracy. To assess the electromagnetic environment created by fixed RF transmitters, an electromagnetic site survey should be considered.

 The presence of RF emitters (such as RFID, MRI, diathermy, electrocautery, etc.) or electronic article surveillance (EAS) in the use environment may cause electromagnetic disturbance in the DolorClast® Focused Shock Waves device operation. Please remove any potential sources of electromagnetic disturbance from the use environment during operation of the DolorClast® Focused Shock Waves device.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT – GUIDANCE
RF conducted IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz outside ISM bands and radio amateur band *	6 Vrms 150 kHz to 80 MHz outside ISM bands and radio amateur band *	Should the measured field strength in the location in which the device is used exceed the applicable RF compliance level above, and the device present abnormal performances, additional measures such as re-orienting or relocating the DolorClast® Focused Shock Waves system may be necessary.
	6 Vrms 150 kHz to 80 MHz in ISM bands and radio amateur band*	6 Vrms 150 kHz to 80 MHz in ISM bands and radio amateur band*	
Proximity fields from RF wireless commu- nication equipment IEC 61000-4-3	27 V/m 380-390 MHz 50% PM 18 Hz 28 V/m 430-470 MHz FM ±5 kHz deviation, 1 kHz sine 9 V/m 704-787 MHz 50% PM 217 Hz 28 V/m 800-960 MHz 50% PM 18 Hz 28 V/m 1700-1990 MHz 50% PM 217 Hz 28 V/m 2400-2570 MHz 50% PM 217 Hz 9 V/m 5100-5800 MHz 50% PM 217 Hz	27 V/m 380-390 MHz 50% PM 18 Hz 28 V/m 430-470 MHz FM ±5 kHz deviation, 1 kHz sine 9 V/m 704-787 MHz 50% PM 217 Hz 28 V/m 800-960 MHz 50% PM 18 Hz 28 V/m 1700-1990 MHz 50% PM 217 Hz 28 V/m 2400-2570 MHz 50% PM 217 Hz 9 V/m 5100-5800 MHz 50% PM 217 Hz	RF wireless equipment maximum output power and separation distance tested (at 30 cm): TETRA 400: Max. 1.8 W GMRS 460, FRS 460: Max. 2 W LTE Band 13, 17: Max. 0.2 W GSM 800/900: Max. 2 W TETRA 800: Max. 2 W iDEN 820: Max. 2 W CDMA 850: Max. 2 W LTE Band 5: Max. 2 W GSM 1800/1900: Max. 2 W CDMA 1900: Max. 2 W DECT: Max. 2 W LTE Band 1,3,4,25: Max. 2 W UMTS: Max. 2 W Bluetooth: Max. 2W WLAN 802.11b/g/n: Max. 2 W RFID 2450: Max. 2 W LTE Band 7: Max. 2 W WLAN 802.11 a/n: Max. 0.2 W Interference may occur in the vicinity of equipment marked with the following symbol: 

*ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 - 6.795 MHz, 13.553 - 13.567 MHz, 26.957 - 27.283 MHz and 40.66 - 40.7 MHz. Amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz - 2 MHz, 3.5 - 4.0 MHz, 5.3 - 5.4 MHz, 7 - 7.3 MHz, 10.1 - 10.15 MHz, 14 - 14.2 MHz, 18.07 - 18.17 MHz, 21.0 - 21.4 MHz, 24.89 - 24.99 MHz, 28.0 - 29.7 MHz and 50.0 - 54.0 MHz.

If the field strength measured in the location in which the DolorClast® Focused Shock Waves system is used exceeds the applicable RF compliance level above, the DolorClast® Focused Shock Waves system should be observed to verify normal operation. If abnormal performances are observed, additional measures may be necessary, such as reorienting or relocating the DolorClast® Focused Shock Waves system

12. ELECTROMAGNETIC COMPATIBILITY (CONT.)

Compatible Cables and Accessories

⚠ The use of accessories and cables other than those specified or sold by EMS as replacement parts may result in increased emissions or decreased immunity of this product.

CABLES AND ACCESSORIES	MAXIMUM LENGTH	COMPATIBLE WITH
Handpiece	< 2.9 m	CISPR 11 Class B / Group 1: RF electromagnetic disturbance
Mains power cable	> 2.9 m	IEC 61000-4-2 Electrostatic discharge (ESD) IEC 61000-4-3 Electromagnetic fields radiated by radio-frequencies IEC 61000-4-4 Electric fast transient (EFT) / burst IEC 61000-4-5 Surge IEC 61000-4-6 Disturbances induced by radio-frequency fields IEC 61000-4-8 Power frequency magnetic field (50/60 Hz) IEC 61000-4-11 Voltage dips, short interruptions and voltage variations

Essential Performance

The DolorClast® Focused Shock Waves system does not present life-sustaining functions, nor does it offer diagnoses of life-supporting functions.

The following functions are presented:

- Continuous shock wave operation
- No changes of mode
- No changes in setpoint values

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
EMS Electro Medical Systems S.A. | EMS worldwide offices (medical) | Affiliates

SWITZERLAND

EMS Electro Medical Systems S.A.
Chemin de la Vuarpillière, 31
1260 Nyon, Suisse
Tel. +41 22 99 44 700
Fax +41 22 99 44 701
e-mail: welcome@ems-ch.com

Manufacturer

EMS Electro Medical Systems S.A.
Chemin de la Vuarpillière 31
1260 Nyon
SWITZERLAND

EC Importer
EC Authorized representative
EMS France SARL
32 Rte de Pontarlier
39460 Foncine-Le-Haut
FRANCE

FRANCE

EMS France Sarl
32, Route de Pontarlier
F-39460 Foncine-le-Haut
Tél. +33 1 84 51 90 01
Fax +33 1 84 51 94 00
e-mail: info@ems-france.fr

ITALY

EMS Italia S.r.l
Via Pisacane 7/b
I-20016 Pero (MI)
Tel. +39 02 3453 8075
Fax +39 02 3453 1724
e-mail: medical@ems-italia.it

GERMANY

EMS Electro Medical Systems GmbH
Stahlgruberring 12
D-81829 München
Tel. +49 89 42 71 61 0
Fax +49 89 42 71 61 60
e-mail: info@ems-ch.de

CHINA

医迈斯电子医疗系统贸易(上海)有限公司
EMS Trading (Shanghai) Co., LTD
24A, Jin Sui Mansion,
N°379 Pudong Nan Rd.
Shanghai, China 200120
Phone: +86 21 3363 2323
Fax +86 21 3363 2152
e-mail: emschina@ems-ch.com

JAPAN

EMS Japan branch Office
Eminent Ningyoucho Bld 6F.,
9-18, Kobuna-cho, Nihombashi,
Chuo-ku, Tokyo, 103-0024, Japan
Tel number +81 3 6810 9351
Fax number +81 3 6810 9352
e-mail: emsjapan@ems-ch.com

JORDAN

EMS Electro Medical Systems Jordan
No.2 - second floor
195 Arar (Wadi Saqra) str.
JO - Amman

SPAIN

EMS Electro Medical Systems España SLU
c/Tomás Bretón, 50-52, 2ª planta
E-28045 Madrid
Tlf. +34 91 528 99 89
e-mail: info@ems-espana.com

USA / CANADA

Electro Medical Systems L.L.C.
2150E Continental Boulevard
Southlake, Texas, 76092, USA
Tel. +1 972 690 8382
Fax +1 972 690 8981
e-mail: emsrepairs@ems-na.com

SWEDEN

EMS Nordic AB
Hammarby Fabriksväg 61
SE - 120 30 Stockholm
Tel. +46 8 899 102
e-mail: info@ems-nordic.se

UNITED KINGDOM

EMS UK Ltd Suite 2
Xenus House, Sandpiper Court
Phoenix Park - Eaton Socon
UK - St Neots, Cambridgeshire
PE19 8EP
Tel. +44 01480 587260
e-mail: info@ems-unitedkingdom.com

AUSTRALIA

EMS Oceania Pty Ltd
The Woolstores
Suite 3, Shed 73
4E Huntley Street
AU - Alexandria, NSW 2015
Tel. +61 (02) 9313 4392
e-mail: info@ems-oceania.com

POLAND

EMS Poland SP. Z O.O.
Al. Roździeńskiego 188H
PL - 40-203 Katowice
Tel. +48 32 4937060
e-mail: info@ems-poland.com

SOUTH KOREA

EMS Electro Medical Systems Korea Ltd
13, Beobwon-ro 8-gil, Songpa-gu
Tel. +82 70 7807 2875
e-mail: emskorea@ems-ch.com

© Copyright EMS S.A. FB-1012/11 rev.B ed.2024/09

Signature Manifest

Document Number: FB-1012/11

Title: DolorClast ® Focus Shock Waves IFU

Effective Date: 14 Oct 2024

Revision: B

All dates and times are in UTC+01:00.

DDM - Release IFU DolorClast Focused Shock Waves (FB-1012) MDR

Refus de collaboration

Name/Signature	Title	Date	Meaning/Reason
Ferran Hernandez (FERNANDEZ)		10 Oct 2024, 02:40:43 PM	Complete

Author and Reviewers

Name/Signature	Title	Date	Meaning/Reason
Ferran Hernandez (FERNANDEZ)		10 Oct 2024, 02:41:03 PM	Approved

Final Reviewer

Name/Signature	Title	Date	Meaning/Reason
Suzanne Fassio Hardy (SFA)	Regulatory Affairs Group Leader	10 Oct 2024, 04:02:29 PM	Approved

Approval

Name/Signature	Title	Date	Meaning/Reason
Timothée Deblock (TDE)	Head of Quality	14 Oct 2024, 02:51:04 PM	Approved

Update ERP

Name/Signature	Title	Date	Meaning/Reason
Marie-Christine Lamy-Biehler (MCL)		14 Oct 2024, 04:20:25 PM	Approved

Notification

Name/Signature	Title	Date	Meaning/Reason
Generic user Notificator (NOTIFICATOR)		14 Oct 2024, 04:20:26 PM	Email Sent
Ferran Hernandez (FERNANDEZ)		14 Oct 2024, 04:20:26 PM	Email Sent